

Silicone Soft plaster: Study



Study title and source:

Determination of the Adhesive Properties of a First Aid Dressing with Silicone Adhesive in Comparison to an Acrylate-Based Adhesive

*Kuhlmann M. et al., EWMA (European Wound Management Association) Conference, Amsterdam, Netherlands, 2017)

Authors:

Kuhlmann M.¹, Klensang K.², Kristof S.², Buck B.², Krause-Kyora F.¹, Wolber R.¹, Grönniger E.²

¹ Beiersdorf Health Care, Beiersdorf AG, Hamburg, Germany, ² Beiersdorf AG, Research & Development, Hamburg, Germany

Study objective:

To determine the gentle removal of a first aid dressing with silicone adhesive (Hansaplast/ Elastoplast Silicone Soft) in comparison to a traditional acrylate-based adhesive

Test dressing/sheet:

- Hansaplast/Elastoplast Silicone Soft
- Standard first aid dressing with acrylate-based adhesive

Study design:

Single-centre, randomised, controlled, blinded investigation

Participants:


30 volunteers, aged 25-63 with intact, undamaged skin

The plaster can be removed completely painless

The plaster is painless to remove compared with my standard wound plaster


For which product did you experience less pain?

Which product was easier/ more convenient/less painful to remove?




■ Silicone adhesive ■ Acrylate-based adhesive

Figure: (Modified illustration of Fig.7*): Self assessment of product performance and preference, pain-free removal



Acrylate-based adhesive: 82% coverage



Silicone adhesive**: 1.6% coverage

Figure: (Modified illustration of Fig.6*): Microscopy (example), stained corneocytes on adhesive matrix after removal

** Hansaplast/Elastoplast Silicone Soft

Results:

Minimal epidermal skin stripping with a significantly lower coverage of corneocytes on the adhesive matrix of the silicone adhesive compared to the acrylate-based adhesive could be shown. Measurements of transepidermal water loss (TEWL) and skin colour demonstrated a significant lower impairment of skin barrier function and skin irritation for the silicone adhesive. Pain-free removal of the silicone adhesive** was confirmed by the test persons.

Conclusion:

The results demonstrated the gentle and skin-friendly removal of the silicone adhesive** with minimal skin damage and minimal stripping of epidermal layers and without causing pain. For self-treatment of minor wounds, first aid dressings based on silicone adhesives should be taken into consideration, particularly when fragile, sensitive or easily irritated skin is affected.